




Figure 2: BovineSNP50 BeadChip Probe Spacing



The BovineSNP50 BeadChip provides uniform coverage across the entire genome. Fewer than 5% of all gaps are larger than 100kb.

provides an average probe spacing of 49.4kb and a median spacing of 36.9kb (Figure 2). Current research on linkage disequilibrium (LD) in multiple breeds of cattle suggests haplotype blocks of approximately 70kb on average<sup>9</sup>, indicating that the resolution offered by the BovineSNP50 chip is well within the resolution of LD in cattle. More than 54,060 SNP probes on the BovineSNP50 BeadChip map to UMD 3.0, the most current bovine reference genome assembly.

### High-Quality Data

The 54,609 SNP assays on the BovineSNP50 BeadChip were subjected to rigorous functional testing to ensure strong performance using Illumina’s Infinium HD Assay. Whole-genome association studies are successful, in part, due to high call rates and accurately called genotypes. Since complex traits often have relatively small gene effects, potential associations may be missed if the assayed SNP in linkage disequilibrium with the SNP of interest has a low call rate or incorrect genotype call. Illumina ensures that every BovineSNP50 BeadChip offers > 99% call rate.

Table 2 shows the outstanding results from Illumina’s internal validation of the BovineSNP50 BeadChip content using samples provided by the Bovine HapMap Consortium. Product developers identified and retained 4,290 loci that appear to have an adjacent or underlying deletion among the breeds sampled. Although these loci yield lower call rates when compared to the other loci on the panel, they were retained because they may provide biologically relevant information for traits of interest and for future improvements in the genome build.

The performance and content validation results clearly show the reliability and outstanding data quality the BovineSNP50 BeadChip delivers. With this BeadChip, researchers can predict the genetic merit for phenotypes of interest and investigate the genetic basis of variation among a multitude of cattle breed types.

### Solutions for Genotyping

As with all of Illumina’s standard products, an optional Laboratory Information Management System (LIMS) and robotic automation are available to accurately and efficiently track samples throughout analysis. Researchers can genotype their own samples using the Illumina GenomeStudio® DNA analysis software or Illumina’s FastTrack Genotyping Service.

Table 1: BeadChip Content Sources

Source	BovineSNP50 v1 Probes	BovineSNP50 v2 Probes
Novel SNPs derived from sequencing common cattle breeds using Illumina’s Genome Analyzer	23,840	24,181
Bovine HapMap Data Set	12,298	12,342
Btau assembly SNPs	9,361	9,404
Whole-Genome Shotgun Reads*	5,808	6,038
Holstein BAC Sequence Data	1,409	1,411
Parentage**	116	120
Other†	1,169	1,113
<b>Total</b>	<b>54,001</b>	<b>54,609</b>

\* Reads derived from six cows/breeds (Norwegian Red, Hostein, Brahman, Angus, Jersey and Limousin) compared against Btau2.0

\*\* Highly curated parentage markers that include 30 duplicates

† Includes common SNPs validated by the Institute for Food and Agricultural Sciences Alberta and INRA, the French International Institute of Agriculture

## Product Summary

Developed in collaboration with leading bovine researchers, the BovineSNP50 BeadChip features more than 54,000 evenly spaced SNP probes that span the bovine genome. This 24-sample BeadChip presents a high-throughput, cost-effective solution for whole-genome studies in beef and dairy cattle.

**Table 2: Product Performance and Specifications**

Parameter	Results	Product Specification
Average Call Rate*	99.7%	> 99%
Reproducibility	100%	> 99.9%
Mendelian Inconsistencies**	0.06%	< 0.1%

\* Based on 473 individuals from 19 major cattle breeds, 31 trios, and 1 replicate

\*\* Based on 4 Holstein trios

**Table 3: BeadChip Content Validation**

Breed	Samples	Polymorphic Loci*	Mean MAF	Median MAF**
Angus	24	41,073	0.21	0.21
Beefmaster	23	43,114	0.22	0.22
Bos indicus Gir	21	23,567	0.11	0.02
Bos indicus Nelore	19	25,492	0.11	0.03
Brahman	22	29,444	0.13	0.07
Brown Swiss	21	35,971	0.19	0.17
Charolais	19	43,723	0.22	0.21
Guernsey	21	36,748	0.19	0.17
Hereford	24	42,132	0.22	0.23
Holstein	49	42,849	0.22	0.22
Jersey	23	35,346	0.18	0.15
Limousin	39	42,617	0.22	0.21
N'Dama	23	28,869	0.14	0.07
Norwegian Red	17	42,055	0.22	0.21
Piedmontese	21	41,912	0.22	0.21
Red Angus	10	42,388	0.21	0.20
Romagnola	21	38,524	0.20	0.19
Santa Gertrudis	21	41,783	0.22	0.21
Sheko	16	35,084	0.17	0.13
<b>Overall</b>	<b>434</b>	<b>47,168</b>	<b>0.24</b>	<b>0.25</b>

\* MAF > 0.05

\*\* Across all 54,609 loci

## References

1. ftp://ftp.hgsc.bcm.tmc.edu/pub/data/Btaurus/fasta
2. www.bovinehapmap.org
3. Gunderson KL, Steemers FJ, Lee G, Mendoza LG, Chee MS (2005) A genome-wide scalable SNP genotyping assay using microarray technology. *Genet* 37(5): 549-554.
4. Steemers FJ, Weihua Chang W, Lee G, Barker DL, Shen R, et al. (2006) Whole-genome genotyping with the single-base extension assay. *Nat Methods* 3(1): 31-33.
5. ftp://ftp.hgsc.bcm.tmc.edu/pub/data/Btaurus/snp/Btau20070913/README
6. Heaton MP, Keen JE, Clawson ML, Harhay GP, Bauer N, et al. (2005) Use of bovine single nucleotide polymorphism markers to verify sample tracking in beef processing. *J Am Vet Med Assoc* 226(8):1311-1314.
7. Matukumalli LK, Taylor JF, and Van Tassell CP. Personal communication.
8. Van Tassell CP, Smith TPL, Matukumalli LK, Taylor JF, Schnabel, RD, et al. (2008) Simultaneous SNP discovery and allele frequency estimation by high-throughput sequencing of reduced representation libraries. *Nat Meth* (accepted).
9. Khatkar MS, Zenger KR, Hobbs M, Hawken RJ, Cavanagh JAL, et al. (2007) A Primary Assembly of a Bovine Haplotype Block Map Based on a 15,036-Single-Nucleotide Polymorphism Panel Genotyped in Holstein-Friesian Cattle. *Genetics* 176(2):763-772.

Ordering Information

Table with 3 columns: Catalog No., Product, and Description. It lists three kits: WG-450-2001 (2 BeadChips), WG-450-2002 (12 BeadChips), and WG-450-2003 (48 BeadChips).

Each BovineSNP50 v2 DNA Analysis BeadChip can process 24 samples and analyze approximately 54,000 loci.

FOR RESEARCH USE ONLY

© 2010 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CPro, GenomeStudio, Genetic Energy, HiSeq, and HiScan are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Pub. No. 370-2007-029 Current as of 03 May 2010



AGAATGATAACAGTAAACACTTCTGTTAACCTTAAGATTACTTGATCCAAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTCTTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC...